



Global Prostate
Solutions

Free Hand Transperineal Needle Guide



Global Prostate
Solutions

BIOPSY KIT

WWW.globalprostatesolutions.com

Free Hand Transperinal Needle Guide


This device is used to direct a needle or instrument into its proper course when performing a clinical and or surgical procedure with the use of the ultrasound probe . This is a single-use device.

Key Features

- ▶ Stable needle alignment
- ▶ Minimises deviation for reliable results
- ▶ Compatible with multiple needle sizes

Specifications

- ▶ Coaxial passes through Full Hole
- ▶ Material: ABS
- ▶ Function: Needle alignment for biopsies



Our Global Presence



With a growing footprint across Europe, Asia, the Middle East, and Africa, our global presence enables us to deliver advanced prostate biopsy technology and support to diverse healthcare systems. We ensure access to reliable, cost-effective prostate biopsy solutions backed by expert training.



 +91 85473 02216, 80738 51840

 contactus@globalprostatesolutions.com

 www.globalprostatesolutions.com